

Customizable configurations. Superior signal quality.

Keypoint Focus keeps patients your priority in EMG/NCS/EP testing

Building on nearly 60 years of EMG system design expertise, Keypoint Focus is a highly flexible and easily scalable EMG solution renowned for its superior signal quality. Accommodating a wide range of amplifiers and other hardware options, Keypoint Focus offers intuitive operation and precise data acquisition, giving your team greater confidence in the care you deliver.

- High CMRR and signal-to-noise ratio for consistent recordings
- Cart-based or portable configurations with 3- to 8-channel amplifiers
- Fast, automatic report generation



Outstanding signal quality

Keypoint's industry-leading amplifiers and stimulators feature outstanding signal quality and reliability. Choose 3-, 4-, 6- or 8-channel system with dedicated inputs for EMG, NCS and EP recordings using either needle electrodes or surface electrodes.

- High CMRR and low input noise provides consistent, high quality recordings
- Software controlled interconnection of reference inputs
- Electrode impedance measurement with LED feedback

Flexible hardware configurations

Keypoint Focus lets you stay closer to patients – regardless of your clinical environment. Choose from a convenient and highly mobile laptop setup, a cart-based system or dedicated workstation, as well as your choice of optional hardware for intuitive operation – including advanced stimulators.

- Available in cart-based or portable configurations
- Advanced stimulus probe with multi-controls for stimulus intensity, stimulus activation, single/recurrent stimulation, reverse anode and cathode and software function shortcuts
- Optional footswitch controls for hands-free operation of system and stimulator
- Flexible amplifier/stimulator arm for close patient connection – easily moved without tools for placement on either side of the system
- InVisus Pro neuromuscular ultrasound with optional mounting hardware for portable use and workstations




Versatile EMG/NCS/EP software

Delivering the best patient care depends on a system that offers flexibility to meet your team's unique clinical needs. Now powered by Keypoint, Natus Elite is the advanced EMG software users trust to assess patients easily and accurately – and quickly create comprehensive reports.

- Industry-leading algorithms minimize study times and improve editing capabilities for simplified review and report generation
- Intuitive interface flexible configuration options
- New testing features, including MUNIX, CMAP scan and ECG-SSEPs enhance patient assessment
- Seamless NMUS integration to complement EMG studies with top-quality imaging
- Automatic temperature measurement for Nerve Conduction studies

Adaptable networking

Keypoint's network capabilities were developed to support a wide range of installation sites while focusing on security and reliability. This makes Focus easily adaptable to small clinics with no professional IT support, as well as large hospital installations with system access controlled by IT using Active Directory Services.



To learn more about Natus products,
contact your local distributor or sales representative.
International Customers Call: **+1-608-829-8500**

Making sense of the body's signals

©2025 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. **046861 RevB**

natus.

Natus Medical Incorporated

natus.com