



Robust, reliable, and advanced EMG

Delivering superior EMG/NCS/EP signal quality to support your diagnostic accuracy

Designed for premium performance, Keypoint G4 is a comprehensive EMG/ NCS/EP solution with renowned signal quality. Accommodating a wide range of amplifiers and other hardware options, Keypoint G4 offers intuitive operation and precise data acquisition, enabling your care team to confidently perform patient assessment and research.

- Convenient, all-in-one workstation with integrated loudspeaker and dedicated controls
- Industry-leading signal quality with high CMRR signal-to-noise ratio
- Scalable amplifiers to meet your unique clinical needs



Outstanding signal quality

Keypoint's industry-leading amplifiers and stimulators feature outstanding signal quality and reliability. Choose 3-, 4-, 6- or 8-channel system with dedicated inputs for EMG, NCS and EP recordings using either needle electrodes or surface electrodes.

- High CMRR and low input noise provides consistent, high quality recordings
- Software controlled interconnection of reference inputs
- Superior automated marker placement
- Electrode impedance measurement with LED feedback

Convenient and comprehensive workstation

Maximize your EMG space with this all-in-one, ergonomically designed workstation, with many adjustable features to keep comfort a priority during every examination.

- Central cart console contains integrated loudspeaker for real EMG sound, high performance ultra-small-form-factor PC and all essential cables
- Flexible amplifier/stimulator arm for close patient connection – easily moved without tools for placement on either side of the system
- Advanced stimulus probe with multi-controls for stimulus intensity, stimulus activation, single/recurrent stimulation, reverse anode and cathode and software function shortcuts
- Optional footswitch controls for hands-free operation of system and stimulator
- InVisus Pro neuromuscular ultrasound with optional mounting hardware




Versatile EMG/NCS/EP software

Delivering the best patient care depends on a system that offers flexibility to meet your team's unique clinical needs. Now powered by Keypoint, Natus Elite is the advanced EMG software users trust to assess patients easily and accurately – and quickly create comprehensive reports.

- Industry-leading algorithms minimize study times and improve editing capabilities for simplified review and report generation
- Intuitive interface flexible configuration options
- New testing features, including MUNIX, CMAP scan and ECG-SSEPs enhance patient assessment
- Seamless NMUS integration to complement EMG studies with top-quality imaging
- Automatic temperature measurement for Nerve Conduction studies

Adaptable networking

Keypoint's network capabilities were developed to support a wide range of installation sites while focusing on security and reliability. This makes the G4 easily adaptable to small clinics with no professional IT support, as well as large research centers and hospital installations with system access controlled by IT using Active Directory Services.



To learn more about Natus products,
contact your local distributor or sales representative.
International Customers Call: **+1-608-829-8500**

Making sense of the body's signals

©2025 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. **046860 RevB**

natus.

Natus Medical Incorporated

natus.com