

THE BASIC SOLUTION FOR REPETITIVE TRANSCRANIAL MAGNETIC STIMULATION

MagPro R20



- ✓ The optimal choice when running a limited number of TMS sessions per day
- ✓ CE approved for depression treatment*
- ✓ Quick and simple protocol set-up
- ✓ Mobile

NEED MORE OPTIONS?

The MagPro R20 Wi-Fi package contains a secure Wi-Fi connection and a configurator, providing you with more options:

- Protocol set up on an external computer
- Wireless upload of the protocol and settings on the MagPro R20
- Print documentation of performed sessions with your logo, header and patient details



Technical data



Output power roll-off	
Max. output at 5 pps	100 %
Max. output at 10 pps	76 %
Max. output at 15 pps	61 %
Max. output at 20 pps	54 %
Mechanical Data	
HxWxD	150 x 390 x 440mm 5.9 x 15.3 x 17.3 in
Weight	20kg (44 lbs)
Technical specifications	
Maximum Repetition Rate	20 pps
Pulse mode	Standard
Waveform	Biphasic
Current direction	Normal
Ordering Number	
9016E0861	



Your local contact:



www.magventure.com

***MagVenture TMS Therapy® is CE approved for:**
Treatment of Major Depressive Disorder in adult patients who have failed to achieve satisfactory improvement from two prior antidepressant medications, at or above the minimal effective dose and duration in the current episode.

For subjects responding to the treatment of MDD,

who have symptoms of comorbid anxiety, the clinical benefit can additionally include an improvement of anxiety symptoms.

Treatment of psychoactive substance use disorder (PSUD) in adult patients.

An adjunct for the treatment of adult patients suffering from Obsessive-Compulsive Disorder (OCD).