

THE ALL-ROUND SOLUTION FOR TRANSCRANIAL MAGNETIC STIMULATION

MagPro R30



- ✓ Quick protocol set-up, storage and recall of predefined set-ups
- ✓ Easily connected to external equipment
- ✓ Data capture and reporting functionality
- ✓ Compatible with all coils
- ✓ Enhanced stimulation capabilities with MagOption
- ✓ CE approved for depression treatment*

NEED MORE OPTIONS?

The MEP Monitor is mounted on the back of MagPro:

- One channel EMG amplifier
- Optimized for Motor Evoked Potentials recordings
- Post processing of MEP Recordings
- Import to spreadsheet program



Technical data



	MagPro R30	MagOption
Mechanical Data		
HxWxD	210 x 530 x 400mm 8.3 x 20.9 x 15.7 in	130 x 530 x 400mm 5.1 x 20.9 x 15.7 in
Weight	36 kg (79 lbs)	28 kg (62 lbs)
Technical specifications		
Maximum Repetition Rate	30 pps, 60 pps, 80 pps	30 pps
Pulse mode	Standard	Standard, Dual and Twin
Waveform	Biphasic, Biphasic Burst (Theta Burst)	Monophasic + Biphasic
Current direction	Normal	Normal
Ordering Number		
	9016E0721	9016E0741



Your local contact:



www.magventure.com

***MagVenture TMS Therapy® is CE approved for:**
Treatment of Major Depressive Disorder in adult patients who have failed to achieve satisfactory improvement from two prior antidepressant medications, at or above the minimal effective dose and duration in the current episode.

For subjects responding to the treatment of MDD,

who have symptoms of comorbid anxiety, the clinical benefit can additionally include an improvement of anxiety symptoms.

Treatment of psychoactive substance use disorder (PSUD) in adult patients.

An adjunct for the treatment of adult patients suffering from Obsessive-Compulsive Disorder (OCD).